


SPEED CAN KILL, BUT IT CAN ALSO SAVE LIVES.

“ In attendance at a motorcycle race, a rider was thrown off his vehicle at high speed and suffered an open complicated fracture with femoral laceration, which was extremely challenging to treat. Due to the nature of the wound, a tourniquet wasn't an option, so as soon as I could get to the casualty, I used Celox. This stopped the bleeding fast and kept the wound sealed, despite a bumpy ride by air ambulance. Without Celox to hand, the story could have ended very differently. ”

FOR MORE INFORMATION VISIT CELOXMEDICAL.COM



*Celox Rapid is indicated to stop life threatening bleeding with a compression time of 60 seconds or until bleeding stops. Hoggarth A, et al. Testing a new gauze hemostat with reduced treatment time. Poster presentation at ATAOC, Ft Lauderdale, August 2011.

The versatile, fast-acting solution to bleeding injuries designed for EMS providers.





セロックスEMSガーゼ 4" x 4"

セロックスEMSガーゼ 8" x 8"

品番: A04053MT
止血ガーゼ 2枚入り

品番: A04054MT
止血ガーゼ 4枚入り



一般医療機器 届出番号:
14B3X10037000005

使用期限: 製造後5年

Specification

特徴

- ・死因の多くを占める大量失血に効果的な止血剤セロックスを含んだガーゼ
- ・エビの殻を精製して抽出したキトサンを使用止血剤としての効果を認められた医療機器 (FDA, CE認可) アレルギー要素は精製時に除去
- ・セロックスEMS止血ガーゼを傷口に適用し、圧迫します。キトサンが血液と反応し、粘着性ゲルを生成。傷口を覆い出血を止めます
- ・傷口に継続して圧迫をかけるため、ガーゼやバンドエーजで傷口を併用すると効果的です
- ・体の止血機構に非依存、非発熱です

サイズ

- セロックスEMSガーゼ 4" x 4"
展開時: 10.2cm x 20.4cm (ガーゼ2枚入り)
- セロックスEMSガーゼ 8" x 8"
展開時: 20.4cm x 20.4cm (ガーゼ4枚入り)
- ガーゼ (1枚あたり) : 10.2cm x 10.2cm

日本総輸入販売元



株式会社オーストリッチインターナショナル
〒222-0033 横浜市港北区新横浜1-14-20
光正第2ビル301
TEL: 045-470-9041 (代) FAX: 045-470-9043